

Inspection Form

Stryker Power-PRO Series Stretcher

Authorized **stryker** Distributor



Serial Number: _____

Vehicle Number: _____

<p>Reason For Service</p> <p><input type="radio"/> New / Initial Inspection</p> <p><input type="radio"/> Routine Maintenance</p> <p><input checked="" type="radio"/> Refurbishment</p> <p><input type="radio"/> Other</p> <p>Cylinder</p> <p><input checked="" type="radio"/> All fasteners are secure.</p> <p><input checked="" type="radio"/> Inspect for and verify that there are no hydraulic fluid (red) leaks.</p> <p><input checked="" type="radio"/> Extend the cylinder rod completely and wipe down with soft cloth and household cleaner.</p> <p><input checked="" type="radio"/> Inspect the fittings and tighten as necessary.</p> <p>Hydraulics</p> <p><input checked="" type="radio"/> Inspect motor mount and verify that all fasteners are secure.</p> <p><input checked="" type="radio"/> Verify that there are no hydraulic fluid leaks.</p> <p><input checked="" type="radio"/> Inspect the reservoir and verify that there are no leaks.</p>	<p>Cables/Wires</p> <p><input checked="" type="radio"/> Verify there is no damage or pinching of wiring harness, cable, or lines.</p> <p><input checked="" type="radio"/> Verify there are no damaged connectors.</p> <p><input checked="" type="radio"/> Check grounds and connections; verify there are no hanging wires.</p> <p>Manual Backup Release Handle</p> <p><input checked="" type="radio"/> Verify the base extends/ retracts smoothly when the manual back-up release handle is engaged.</p> <p><input checked="" type="radio"/> With 100 lbs or more on the cot, verify the cot does not lower when the manual backup release handle is pulled.</p> <p>Litter</p> <p><input checked="" type="radio"/> Verify all fasteners are secure.</p> <p><input checked="" type="radio"/> Verify backrest operates properly.</p> <p><input checked="" type="radio"/> Adjust pneumatic cylinder for full range motion if required.</p> <p>Base</p> <p><input checked="" type="radio"/> Verify all fasteners are secure.</p> <p><input checked="" type="radio"/> Inspect the cot frame/base.</p>	<p>X-Frame</p> <p><input checked="" type="radio"/> Verify all fasteners are secure.</p> <p>Head Section</p> <p><input checked="" type="radio"/> Verify all fasteners are secure.</p> <p><input checked="" type="radio"/> Verify the head section extends and locks properly.</p> <p><input checked="" type="radio"/> Inspect the SMRT battery housing and pole mount area for cracks or damage.</p> <p>Wheels</p> <p><input checked="" type="radio"/> Inspect all wheels for wear, rolling and spinning properly.</p> <p>Accessories</p> <p><input checked="" type="radio"/> Verify the I.V. pole operates properly.</p> <p>Labels</p> <p><input checked="" type="radio"/> Good</p>	<p>Inspector Evaluation</p> <p><input checked="" type="radio"/> REFURBISHED – Any faulty or broken parts and components have been replaced with Stryker OEM parts. Equipment is in functional condition and is ready for use. Follow regularly-scheduled maintenance inspections.</p> <p>Acknowledgement of Maintenance, Service, and Repair</p> <p>The routine maintenance, service, or repair was completed and the equipment was evaluated with recommendations noted.</p> <p>Date of Service/Repair: _____</p> <p>Signature of Service/Repair Person: _____</p>
---	---	--	---