

Rowland Incubator Adaptor Cot Compatibility List

The Rowland Incubator Adapter is a device that provides a stable location on a cot to place an incubator deck. This document applies to both the R-105041 and R-105042 kits. The incubator deck must have the post and wheel cup layout of a Ferno #9 cot on its underside. In addition to the Incubator Adapter, certified straps must be used to secure the deck to the cot. Below is a list of popular EMS cots indicating what models are compatible.

Manufacturer	Model	Serial Number	Compatible
Stryker	6500 Power Pro	060939080 and up	Yes
Stryker	6506 Power Pro	All	Yes
Stryker	6080 MX-Pro	All	Yes
Stryker	6082 MX-Pro	All	Yes
Ferno	All	All	No

